

Recommended Shelf Life

ROMIL-UpS™ Ultra Purity Solvents

Code	Description	Storage	Shelf Life	Opened
H033	Acetone UpS ultra gc	RT	2 years	3 months
H054	Acetonitrile UpS ultra gc	RT	2 years	
H050	Acetonitrile UpS ultra lc	RT	2 years	
H204	Dichloromethane UpS stabilised with amylene ultra lc	RT	1 year	3 months
H205	Dichloromethane UpS stabilised with cyclohexene/amylene ultra gc	RT	1 year	3 months
H317	Ethanol absolute UpS ultra lc	RT	2 years	
H347	Ethyl Acetate UpS ultra lc	RT	2 years	
H364	n-Heptane 99% UpS ultra gc	RT	2 years	
H391	n-Hexane 95% UpS ultra lc	RT	2 years	
H394	n-Hexane 99% UpS ultra gc	RT	2 years	
H412	Methanol UpS purge and trap	RT	2 years	
H415	Methanol UpS ultra gc	RT	2 years	
H411	Methanol UpS ultra lc	RT	2 years	
H573	n-Pentane 99% UpS ultra gc	RT	2 years	
H626	Propan-2-ol UpS ultra lc	RT	2 years	
H720	Tetrahydrofuran UpS ultra lc	RT	2 years	3 months
H772	Toluene UpS ultra gc	RT	2 years	
H903	2,2,4-Trimethylpentane UpS ultra gc	RT	2 years	
H949	Water UpS ultra lc	RT	2 years	1 month

Code	Storage Temperature & Conditions
RT	Store at 15-25°C (room temperature)
RD	Store at 15-25°C in dark (room temperature)
FR	Store at 2-8°C (fridge)
FD	Store at 2-8°C in dark (fridge)
FZ	Store at -10 to -25°C (freezer)
DI	Store at -70°C (dry ice)

Shelf Life period after opening must run within the unopened period otherwise it should be reduced accordingly.

For more information see www.romil.com/products/shelf-lives/

Published: 1 March 2016