

## Recommended Shelf Life

### ROMIL-UpS™ Ultra Purity Solvent Mixes

Code	Description	Storage	Shelf Life	Opened
R8445	LC Eluant mix AA2 UpS acetonitrile with acetic acid 0.1% ultra lc	RT	1 year	
R5803	LC Eluant mix AF2 UpS acetonitrile with formic acid 0.1% ultra lc	RT	1 year	
R8384	LC Eluant mix AT4 UpS acetonitrile with trifluoroacetic acid 0.1% ultra lc	RT	1 year	
R5604	LC Eluant mix MA2 UpS methanol with acetic acid 0.1% ultra lc	RT	1 year	
R5729	LC Eluant mix MF4 UpS methanol/water 90:10 / buffer ultra lc	RT	1 year	
R9038	LC Eluant mix MF5 UpS methanol/water 10:90 / buffer ultra lc	RT	1 year	
R6615	LC Eluant mix MF7 UpS methanol with formic acid 0.1% ultra lc	RT	1 year	
R8127	LC Eluant mix MFT1 UpS methanol/water 10:90 / buffer with tfa 0.1% ultra lc	RT	1 year	
R3124	LC Eluant mix MT2 UpS methanol with trifluoroacetic acid 0.1% ultra lc	RT	1 year	
R7416	LC Eluant mix WA4 UpS water with acetic acid 0.1% ultra lc	RT	1 year	3 months
R4905	LC Eluant mix WF2 UpS water with formic acid 0.1% ultra lc	RT	1 year	3 months
R2614	LC Eluant mix WT4 UpS water with trifluoroacetic acid 0.1% ultra lc	RT	1 year	3 months
R2257	LC Wash solution AMPW1 UpS acetonitrile/methanol/propan-2-ol/water 25:25:25:25	RT	1 year	
R1485	LC Wash solution PW1 UpS propan-2-ol/water 50:50	RT	1 year	

Code	Storage Temperature & Conditions
RT	Store at 15-25°C (room temperature)
RD	Store at 15-25°C in dark (room temperature)
FR	Store at 2-8°C (fridge)
FD	Store at 2-8°C in dark (fridge)
FZ	Store at -10 to -25°C (freezer)
DI	Store at -70°C (dry ice)

Shelf Life period after opening must run within the unopened period otherwise it should be reduced accordingly.

For more information see [www.romil.com/products/shelf-lives/](http://www.romil.com/products/shelf-lives/)

Published: 1 March 2016